



# POST-OPERATIVE SHOE

With forefoot discharge for percutaneous surgery

## CHAUSSURE PERCUSOFT



- Immediate functional recovery without the use of canes.
- Surgical results protection even with pin in metatarsal phalangeal joint
- Stability and comfort.

- Upper Shoe : Breathable fabric
- Sizes : 35/36 — 37/38 — 39/40 — 41/42 — 43/44 — 45/46.
- Right and Left shoe available

## TECHNICAL SPECIFICATIONS

- The tipping point position combined with the sole stiffness allow a roll of the foot without metatarso phalangeal flexion therefore reducing constraint on the first ray. The selected materials provide important walking shocks absorption without compromising the needed stiffness.
- The upper shoe dimensions have been studied to adapt perfectly to operated foot shape, taking into account post-operative bandages and post-operative swelling. Thanks to the elastic fabric, the flaps can easily fit on any volume allowing a precise fastening. It results in a slimmer and more elegant footwear.
- Building on our Barouk Shoe experience, the Percusoft is available in right and left side separately. It allows a perfect foot holding and a better walking experience.

## INDICATIONS

- Foot surgery, especially for forefoot percutaneous surgery.
- Each acute state with the need of a rigid sole enabling the walking process without forefoot constraints.

*Romans Industrie*